



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Refrigerated Micro Centrifuge, Digital, 16000 rpm

Product Image



Description

Table top refrigerated micro centrifuge machine. The temp. from room temp. upto minus 8°C is controlled by PID electronic digital temp. controller cum indicator. Microprocessor based Digital Timer of range 0-15 min. with LCD display. Preset speed controls thru microprocessor from 2000 to 16000 r.p.m . Digital display shows Set Time/Run Time /R.P.M. S.S. Bowl Size: 9" dia

Special Features

Microprocessor based Table Top centrifuge.
Table top centrifuge compact and elegant design.
Fitted with vacuum feet.
CFC free refrigerant.
Centrifuge bowl made of stainless steel.

Supply:- complete with following angle rotor heads to accommodate polypropylene micro centrifuge tubes (Ependroff Tubes as well).

Specification

Angle Rotor Head Cap. 12x1.5/2ml
Adopters to accept 1.5ml Eppendrof micro-centrifuge tubes.

Also supplied with head puller, spanner, 12 Ependroff Tubes.

Model

Refrigerated Micro Centrifuge, Digital, 16000 rpm

Catalog No.

JA10220

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.