



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Refrigerated P.R.P. Centrifuge Machine

## Product Image



## Description

---

Model : 60

### Specifications:

- Bench Top universal Refrigerated Centrifuge Machine, all purpose.
- Heavy duty C.R.C. construction finished with powder coating.
- Stainless Steel centrifuge bowl for extra protection and easy cleaning.
- Fitted with well balanced universal motor with having carbon bushes.
- 2 lines of 16 character LCD display of r.p.m., for set and run time.
- Two Programs are feeded in software of machine as per recommended by for P.R.P. Tests.
  - (i) First Program 4 minutes, 3400 r.p.m.
  - (ii) Second Programme 2 minutes, 3500 r.p.m. OR we can programme speed from 500 to 5000 r.p.m. in variation of 100 r.p.m. as per Doctors requirement .
- Temp. is controlled from -10°C upto 37°C in intervals of 1°C.
- P.I.D. Electronic digital temperature controller cum indicator is fitted to control the centrifuge temp.

### Rotors:

It will be supplied with two different rotors details as under:

- i) 8 Holes of 15 ml Swing Out Rotor.
- ii) 4 Holes to hold 4 Dr. P.R.P. Syringe of 20 ml each.

-

Refrigerated P.R.P. Centrifuge Machine

**Catalog No. 100756**

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.