



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Refrigerated P.R.P. Centrifuge Machine

Product Image



Description

Model : 60

Specifications:

- Bench Top universal Refrigerated Centrifuge Machine, all purpose.
- Heavy duty C.R.C. construction finished with powder coating.
- Stainless Steel centrifuge bowl for extra protection and easy cleaning.
- Fitted with well balanced universal motor with having carbon bushes.
- 2 lines of 16 character LCD display of r.p.m., for set and run time.
- Two Programs are feeded in software of machine as per recommended by for P.R.P. Tests.
(i) First Program 4 minutes, 3400 r.p.m.
(ii) Second Programme 2 minutes, 3500 r.p.m. OR we can programme speed from 500 to 5000 r.p.m. in variation of 100 r.p.m. as per Doctors requirement .
- Temp. is controlled from -10°C upto 37°C in intervals of 1°C.
- P.I.D. Electronic digital temperature controller cum indicator is fitted to control the centrifuge temp.

Rotors:

It will be supplied with two different rotors details as under:

- 8 Holes of 15 ml Swing Out Rotor.
- 4 Holes to hold 4 Dr. P.R.P. Syringe of 20 ml each.

-

Refrigerated P.R.P. Centrifuge Machine

Catalog No. 100756

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.