



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Research Stereo Microscope with Turret Mount And Dual Illumination, Monoscope 20x

Product Image



Description

- Body:** Binocular head inclined at 45 degrees, 360 degrees rotatable head with adjustable inter-pupillary distance.
- Base:** Stage with clips and transparent glass plate
- Focussing:** Coarse focussing using diagonally cut rack and pinion arrangement.
- Illumination:** With integral illuminators, two 6V/20W lamps are fitted to provide both transmitted and incident illumination. Built-in transformer inside the base. Lights can be used individually or simultaneously.
- Eye Piece:** Wide Field Eyepiece 10X.
- Objective:** 2 sets of objectives 2X and 4X, selected by rotating nosepiece through 90-degrees.
- Box:** Complete in thermocol box with essential accessories.

Model

Catalog No.

Research Stereo Microscope, with Turret Mount and Dual Illumination, Monoscope 20x

JM003800

Note: Available in plywood box on demand on extra cost.

Note: Objectives set of 3X and 5X are available on demand.

Similar to E8R05716

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.