



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Resistivity Of Semiconductors

## Product Image



## Description

---

Resistivity of semiconductors by four probe method at different Temperatures & Determination of Band gap.

**Features :** The experiment consists of the following accessories for four probe arrangement, oven (upto 200°C), Germanium crystal mounted, Thermometer (0-200°C), four probe setup, constant current power supply, digital electronic millivoltmeter of range 200 MV. Display 3 1/2 digit 7 segment LED.

Catalog No.	Particulars
JE20040	Resistivity Of Semiconductors

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.