



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Resistivity of Semiconductors

Product Image



Description

Resistivity of semiconductors by four probe method at different Temperatures & Determination of Bandgap

Features : The experiment consists of the following accessories for four probe arrangement, oven (upto 200oC), Germanium crystal mounted, Thermometer (0-200oC), four probe setup, constant current power supply, digital electronic millivoltmeter of range 200 mV. Display 3 1/2 digit 7 segment LED.

Catalog No.	Particulars
JE20270	Resistivity of Semiconductors

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.