



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Rubber Teats

Product Image



Description

Teats are used in laboratory with droppers for dispensing drops or qualitative amounts of liquids or chemicals.

- Available in Silicon, Natural & Latex Rubber.

Catalog No.	Type	Rubber Type	Capacity
101360-A	Conical	Silicon Rubber	1ml
101360-B	Conical	Silicon Rubber	2-3ml

101360-C	Conical	Natural Rubber	1ml
101360-D	Conical	Natural Rubber	2-3ml
101360-E	Round	Natural Rubber	1ml
101360-F	Round	Latex Rubber	Small
101360-G	Round	Latex Rubber	Large

- 1ml Teat available in Pack of 100.

- 2-3ml Teat available in Pack of 50.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.