



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Serological Water Bath, Digital Temperature Controller, All Stainless Steel

Product Image



Description

Serological water bath is suitable for clinical procedure, incubation, inactivation & pharmaceutical biomedical procedure. Construction: Double walled chambers made of stainless steel & outer wall made of 304 Stainless Steel finished with epoxy based powder coated paint. The double walled chambers is provided with a drain plug to facilitate easy emptying & cleaning of the inner chamber when required S.S. Lid & perforated removable diffuser are also provided with double walled chamber. This chamber house 1500 watts heating device, which is controlled by Digital temperature, Stainless steel outer body from ambient to 100°C within $\pm 1^{\circ}\text{C}$. Two indicator lights & on/off Rotary Switch

Particulars	Catalog No.
Serological Water Bath, Digital Temperature Controller, All Stainless Steel, 300x250x175mm (2 racks)	JA21390
Serological Water Bath, Digital Temperature Controller, All Stainless Steel, 350x300x175mm (4 racks)	JA21391
Serological Water Bath, Digital Temperature Controller, All Stainless Steel, 450x300x175mm (6 racks)	JA21392

Note: All the above models can be fitted with 1/20 HP motor with stainless steel rod & blades.

Note: Customized Baths (Water/Oil) with specific requirements like size, material, temperature range, etc are available on special demand

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.