



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Serological Water Bath, Digital Temperature Controller

## Product Image



## Description

---

Serological water bath is suitable for clinical procedure, incubation, inactivation & pharmaceutical biomedical procedure. Construction: Double walled chambers made of stainless steel & outer wall made of thick steel finished with epoxy based powder coated paint. The double walled chambers is provided with a drain plug to facilitate easy emptying & cleaning of the inner chamber when required S.S. Lid & perforated removable diffuser are also provided with double walled chamber. This chamber house 1500 watts heating device, which is controlled by digital temperature controller cum indicator from ambient to 100°C within  $\pm 1^\circ\text{C}$ . Two indicator lights & on/off Rotary Switch.

Particulars	Catalog No.
Serological Water Bath, Digital Temperature Controller, 300x250x175mm (2 racks)	JA21380
Serological Water Bath, Digital Temperature Controller, 350x300x175mm (4 racks)	JA21381
Serological Water Bath, Digital Temperature Controller, 450x300x175mm (6 racks)	JA21382

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.