



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Serological Water Bath, Thermostatic

Product Image



Description

Serological water bath is suitable for clinical procedure, incubation, inactivation & pharmaceutical biomedical procedure. Construction: Double walled chambers made of stainless steel & outer wall made of thick steel finished with epoxy based powder coated paint. The double walled chamber is provided with a drain plug to facilitate easy emptying & cleaning of the inner chamber when required S.S. Lid & perforated removable diffuser are also provided with double walled chamber. This chamber house 1500 watts heating device, which is controlled by capillary type thermostat from ambient to 100°C within $\pm 1^\circ\text{C}$. Two indicator lights & on/off Rotary Switch.

Particulars

Catalog No.

Serological Water Bath, Thermostatic, 300x250x175mm (2 racks)

JA21370

Serological Water Bath, Thermostatic, 350x300x175mm (4 racks)

JA21371

Serological Water Bath, Thermostatic, 450x300x175mm (6 racks)

JA21372

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.