



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Shaking Water Bath

Product Image



Description

This water bath is suitable for studies on tissue metabolism, enzymes and protein coagulations or other experiments requiring shaking at constant speed under controlled temperature environment.

Construction: The inner chamber consists of thick stainless steel tank housed in a powder coated outer case filled with high-grade insulation to prevent heat loss. Heaters are controlled from ambient to 100°C by capillary type thermostat. This bath has a shaking trolley, which can be controlled from 40 to 140 strokes per minute. The shaking tray can either hold test tubes or flasks (of 25ml or 50ml or 100 ml capacity) as required. The front panel includes ON/OFF switches for heaters & shaking mechanism. It also includes temperature control knob & speed control knob and provided with L shaped thermometer to indicate the inside temperature. Shaking water bath is also available with Digital temperature controller & P.I.D. controller. Shaking Trays of Stainless Steel are Extra as

required

Particulars	Catalog No.
Shaking Water Bath, 275 x 275 x 75 mm with L-shaped thermometer	JA21460
Shaking Water Bath, 275 x 275 x 150 mm with L-shaped thermometer	JA21461
Shaking Water Bath, 275 x 275 x 75 mm with Digital Temperature Controller	JA21462
Shaking Water Bath, 275 x 275 x 150 mm with Digital Temperature Controller	JA21463
Shaking Water Bath, 275 x 275 x 75 mm with P.I.D. Controller	JA21464
Shaking Water Bath, 275 x 275 x 150 mm with P.I.D. Controller	JA21465
Shaking Water Bath, Shaking Tray of Stainless Steel for 16 nos. 25ml Conical Flasks	JA21466
Shaking Water Bath, Shaking Tray of stainless steel for 12 nos. 50 ml Conical Flasks	JA21467
Shaking Water Bath, Shaking Tray of stainless steel for 9 nos. 100 ml Conical Flasks	JA21468
Shaking Water Bath, Shaking Tray of stainless steel for 5 nos. 250ml Conical Flasks	JA21469
Shaking Water Bath, Shaking Tray of stainless steel for 4 nos. 500ml Conical Flasks	JA21470
Shaking Water Bath, Shaking Tray of stainless steel for 2 nos. 1000ml Conical Flasks	JA21471
Shaking Water Bath, Shaking Tray of stainless steel for 25 nos. 25mm dia Test Tube	JA21472

Note: Digital R.P.M. meter also available.

Note: Customized Baths (Water/Oil) with specific requirements like size, material, temperature range, etc are available on special demand

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.