



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Sietz Filter Holders

Product Image



Description

Sietz filters are to remove submicron. Particles for bacteriological and other laboratory work made of Stainless Steel.

Particulars

Catalog No.

Sietz Filter Holders, Stainless steel, 6 cm. pad, 100 ml. capacity

JA23860

Sietz Filter Holders, Stainless steel, 8 cm. pad, 250 ml. capacity

JA23861

Sietz Filter Holders, Stainless steel, 8 cm. pad, 500 ml. capacity

JA23862

Sietz Filter Holders, Stainless steel, 14 cm. pad, 1000 ml. capacity

JA23863

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.