



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Specimen Cutting Machine

## Product Image

---



## Description

---

Specimen cutting is the first step during sample preparation. It adopts high speed rotating thin grinding wheel to cut specimen. It is equipped with cooling system to prevent the heat damage during cutting.

## Specification

|                                 |                        |
|---------------------------------|------------------------|
| <b>Maximum Sample Selection</b> | 30x30 mm               |
| <b>Grinding Wheel Size</b>      | 300x2x32 mm            |
| <b>Rotate Speed</b>             | 1500rpm / min. approx. |
| <b>Dimensions</b>               | 24x26x30 Inches        |
| <b>Net Weight</b>               | 40 Kg                  |

| <b>Model</b>             | <b>Catalog No.</b> |
|--------------------------|--------------------|
| Specimen Cutting Machine | 7001 CM            |

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.