



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Specimen Cutting Machine

Product Image



Description

Specimen cutting is the first step during sample preparation. It adopts high speed rotating thin grinding wheel to cut specimen. It is equipped with cooling system to prevent the heat damage during cutting.

Specification

Maximum Sample Selection	30x30 mm
Grinding Wheel Size	300x2x32 mm
Rotate Speed	1500rpm / min. approx.
Dimensions	24x26x30 Inches
Net Weight	40 Kg

Model	Catalog No.
Specimen Cutting Machine	7001 CM

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.