



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Specimen Mounting Press

Product Image



Description

1. Specimen Mounting Press is manual mechanic mounting press with 3 different size of mounting mould.
2. Temperature and heating time can be setup and shown manual on the panel.

Specification

Mould Diameter	22mm, 130mm, 145mm
Input Voltage	220V/50Hz, 110V/60Hz.
Power	Max.650 W
Dimensions	340x260x430 mm (Approx.)
Weight	33 Kg (Approx.)

Model	Catalog No.
Specimen Mounting Press	7001 SMP

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.