



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Specimen Mounting Press

## Product Image

---



## Description

---

1. Specimen Mounting Press is manual mechanic mounting press with 3 different size of mounting mould.
2. Temperature and heating time can be setup and shown manual on the panel.

### Specification

<b>Mould Diameter</b>	22mm, 130mm, 145mm
<b>Input Voltage</b>	220V/50Hz, 110V/60Hz.
<b>Power</b>	Max.650 W
<b>Dimensions</b>	340x260x430 mm (Approx.)
<b>Weight</b>	33 Kg (Approx.)

<b>Model</b>	<b>Catalog No.</b>
Specimen Mounting Press	7001 SMP

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.