



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Staining Box

Product Image



Description

Material: Polyrpopylene

Moulded in polypropylene this staining box helps in staining fixing de-staining and handling fragile electrophoresis gels and menbranes.Box as well as cover will resist temperature from -70 ° C to + 70° C.The LDPE plug provides leakproof seal and allows easy liquid drainage.

Catalog No.	Size (mm)
JW10101	125x125x50

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.