



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Storage Boxes, Plastic

## Product Image

---



## Description

---

**Material:** Polypropylene

These storage boxes made of polypropylene prove to be an excellent option for storage of different powdered reagents and other material in the laboratory .Square shape of these boxes help in saving shelf space .Available in three different sizes these boxes have polypropylene caps are autoclavable provided good chemicals resistance and heavy duty services.

Capacity (ml)	Catalog No.
250	JW12401
500	JW12402
1000	JW12403

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.