



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Student Organ Bath (Four Unit)

## Product Image



## Description

### THERMOSTATIC ORGAN BATH:

Organ Bath with fine thermostatic temperature control and complete set of accessories, making a versatile unit for all research work.

Acrylic Bath 14"x 9"x 8" with leak proof joints, rustling in sturdy metal frame, inside the bath is metal sheathed heating element, controlled by means of a adjustable thermostat.

Thermostat sensitivity  $\pm 1^{\circ}\text{C}$

Acrylic junction box is fitted to the bath frame, and carries all the electrical leads and connections. Two red lamps indicate functioning of the mains and stirrer. One green lamp indicates functioning of the thermostat.

Standard Components are:

1. Organ tubes 4
2. Oxygen tubes 4
3. Stainless steel capillary lever 4
4. Warming Quail 4
5. X Block 4
6. Simple Lever 4
7. Hook Grip Rod 4

-

Student Organ Bath (Four Unit)

**Catalog No. 100674**

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.