



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Student Organ Bath (Double Unit)

Product Image



Description

A Perspex Organ Bath with fine thermostatic temperature control, electronic stirrer and complete set of accessories, making it a versatile unit for all research work. Bath 12" x 9" x 8", with leak proof joints, resting in sturdy metal frame; top of the bath carries a 1" x 3/8"

thick band, for fixing accessories, inside the bath is a metal sheathed heating element (easily replaceable) controlled by means of a fine adjustable thermostat with sensitivity $\pm 1^{\circ}\text{C}$. The stirrer unit is completely independent and dips in the bath in one corner, fixing on its special stainless steel upright. A Mild Sheet Junction box is fitted to the bath frame and carries all the electrical leads and connections. Two pilot lamps indicate functioning of the thermostat and stirrer.

Standard components are:

- Double set of glass inner vessels
- Oxygen tubes
- Warming coils
- Two stainless steel capillary levers
- Four stainless steel uprights with boss heads
- X block

-

Student Organ Bath (Double Unit)

Catalog No. 100673

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.