



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Tablet Hardness Tester (Pfizer Type)

Product Image



Description

A handy, tough and direct dial, reading type tablet hardness tester, working as pair of pliers, generates breaking force in single grip. The instrument offers accurate and reproducible result much faster than other conventional hardness tester.

Specifications:

Accuracy: ± 0.2 kg Least

Count: 0.2 kg

Safe Testing Force: 20 kg Maximum

Tablet Diameter: 25 mm

-

Tablet Hardness Tester (Pfizer Type)

Catalog No. 100683

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.