



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Test Tube Baskets

Product Image



Description

Material: Polypropylene

Test Tube baskets are made of polypropylene. These autoclavable unbreakable and non-corrosive baskets are ideal for rinsing and draining glassware. These are also useful for storage of test tubes during sterilization or incubation.

Catalog No.	Size (mm)
JW09201	140x120x110
JW09202	160x160x160

Disclaimer

The Product details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.