



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Test Tube Baskets

## Product Image

---



## Description

---

**Material:** Polypropylene

Test Tube baskets are made of polypropylene. These autoclavable unbreakable and non-corrosive baskets are ideal for rinsing and draining glassware. These are also useful for storage of test tubes during sterilization or incubation.

Catalog No.	Size (mm)
JW09201	140x120x110
JW09202	160x160x160

## Disclaimer

---

The Product details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.