



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Test Tube Holder, Metal Wire

Product Image



Description

Plated steel wire holder with finger grip and wire collar. The jaws accept tubes up to 40mm diameter and their alignment is maintained by a wire collar. Overall length about 16cm

Catalog No.	Particulars
JD20290	Test Tube Holder, Metal Wire

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.