



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Test Tube Screw Cap, Plastic

Product Image



Description

Material: Polypropylene

These polypropylene tubes provide excellent option for RIA, coagulation & bacteriology etc. The cap is screwed on to the tube & hence this tube can be filled up to the brim

Catalog No.	Size (mm)
JW06201	12x75
JW06202	12x100
JW06203	16x100

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.