



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Test Tube Stand (Round), Plastic

Product Image



Description

Material: Polypropylene

A handy little space saver this circular 12 place polypropylene test tube stand has four 25 mm dia holes and eight 19 mm dia hole in its top plates. The base has hemispherical wells and 16 vertical pins for drying tubes surrounded by a trough to catch dripping. This autoclavable rack can be easily disassembled for cleaning purpose.

Catalog No.

JW08501

Item

Test Tube Stand (Round)

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.