



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Test Tube with Stopper

Product Image



Description

Material:-Borosilicate Glass

A type Test Tube with Interchangeable Stopper

B type Test Tube, Graduated with Interchangeable Stopper

A type: Test Tube with Interchangeable Stopper

Catalog No.	Size (mm)
JG20501	12x100
JG20502	15x125
JG20503	18x150
JG20504	22x150
JG20505	25x150
JG20506	25x200
JG20507	32x200
JG20508	38x200

B type: Test Tube, Graduated with Interchangeable Stopper

Catalog No.	Size (mm)	Capacity (ml)
JG20601	12x100	8
JG20602	15x125	10
JG20603	18x150	20
JG20604	22x150	25
JG20605	25x150	50
JG20606	25x200	75
JG20607	32x200	100

All tubes fitted with standard stopper. For custom interchangeable stopper please specify at the time quotation/order

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.