



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Thermometer, Precision

Product Image



Description

Glass capillary bulb tube with printed scale. Length 450mm

Different types available :

Thermometer, Precision, Yellow back, Alcohol

Thermometer, Precision, Yellow back, Mercury

Catalog No.Particulars

JD14740	Thermometer, Precision, Yellow back, Alcohol, -10 to 50 C, 0.2
JD14741	Thermometer, Precision, Yellow back, Alcohol, -10 to 50 C, 0.1
JD14742	Thermometer, Precision, Yellow back, Alcohol, -10 to 110 C, 0.2
JD14743	Thermometer, Precision, Yellow back, Alcohol, -10 to 110 C, 0.1
JD14744	Thermometer, Precision, Yellow back, Alcohol, -50 to + 50 C, 1.0
JD14745	Thermometer, Precision, Yellow back, Alcohol, -50 to + 100 C, 1.0
JD14746	Thermometer, Precision, Yellow back, Alcohol, -100 to + 50 C, 1.0
JD14747	Thermometer, Precision, Yellow back, Alcohol, -100 to + 100 C, 1.0
JD14748	Thermometer, Precision, Yellow back, Mercury, -10 to 50 C, 0.2
JD14749	Thermometer, Precision, Yellow back, Mercury, -10 to 50 C, 0.1
JD14750	Thermometer, Precision, Yellow back, Mercury, -10 to 110 C, 0.2
JD14751	Thermometer, Precision, Yellow back, Mercury, -10 to 110 C, 0.1

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.