




Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Thermometer, Wall, Plastic Body

Product Image



Description

Red spirit filled thermometer, Consist of specially designed glass capillary which facilitates the reading inside column. Having provision for wall mounting. Range: -10 to 50 C

Catalog No.	Particulars
JD14780	Thermometer, Wall, Plastic Body, Small Size
JD14781	Thermometer, Wall, Plastic Body, Big Size

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.