




**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Thermometer, Wall, Plastic Body

## Product Image



## Description

---

Red spirit filled thermometer, Consist of specially designed glass capillary which facilitates the reading inside column. Having provision for wall mounting. Range: -10 to 50 C

Catalog No.	Particulars
JD14780	Thermometer, Wall, Plastic Body, Small Size
JD14781	Thermometer, Wall, Plastic Body, Big Size

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.