



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Thermometer Oral (Prismatic)

Product Image



Description

Thermometer Oral, prismatic type, clinical is a economy model with lens front, dual scale and cylindrical bulb. Available with plastic case.

Range: 35-45 °C & 95-108°F

Accuracy: ± 0.1 °C & ± 0.2 °F

Catalog No.	Particulars
JH10350	Thermometer Oral (Prismatic)

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.