




Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Thermometers

Product Image



Description

General purpose glass thermometers white/yellow glass with screen printing mercury filled packed in cardboard/plastic case total length 300 mm.

General purpose glass thermometers white/yellow glass with screen printing red alcohol filled packed in card board/plastic case total length 300 mm.

Catalog No.	Particulars
--------------------	--------------------

JZ03220	Thermometers, -10 to 1100c Alcohol, each in Card Board case
JZ03221	Thermometers, -10 to 1100c Mercury, each in Card Board Case
JZ03222	Thermometers, -10 to 1100c Alcohol, each in Plastic Case
JZ03223	Thermometers, -10 to 1100c Mercury, each in Plastic Case
JZ03224	Thermometers, 0 to 1500c Mercury, each in Plastic Case
JZ03225	Thermometers, 0 to 2500c Mercury, each in Plastic Case
JZ03226	Thermometers, 0 to 3600c Mercury, each in Plastic Case

A quality product from JAPSON ideal for use in School Laboratories in developing countries of Asia, Africa and South America.

Standard specifications for school laboratory use. Ideal mix of Quality and Economy. Other configurations and superior quality available for specific use or for developed countries.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.