



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Thin Layer Chromatograph (T.L.C. Kit)

Product Image



Description

This unit is based on Dosage German pattern having movable applicator with in-built thickness arrangement between 0 to 2mm. The components are :-

- (a) Spreader (Applicator) made of electroplated brass.
- (b) Perspex Base size 14×23 cm to support 5 glass plates of size 20×5 cm.
- (c) Plate Rack Aluminum for ten 20×20 cm plates.
- (d) Spotting Term plate made of Perspex.
- (e) Developing tank with lid.
- (f) T.L.C. plates set of five 20×20 cm & two 20×5 cm or set of 20×10 cm and two 20×5 cm.
- (g) Micro-pipette.
- (h) Scriber for making lines made of stainless steel.
- (i) Glass sprayer with rubber bellows capacity 100 ml.
- (j) Instruction Manual.
- (b) It is same as above but the applicator is made of anodized aluminum.

ACCESSORIES for T.L.C. Apparatus are also available

-

Thin Layer Chromatograph (T.L.C. Kit)

Catalog No. 100685

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.