



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Tilt Measures

Product Image





Description

Type A Glass with acid proof Nylon stopper (No. 4, 25mm)

Type B Plastic with acid proof Nylon stopper (No. 4, 25mm)

A type: Glass with acid proof Nylon stopper (No. 4, 25mm)

Catalog No.	Capacity (ml)
JG27801	1 ml
JG27802	10 ml

B type: Plastic with acid proof Nylon stopper (No. 4, 25mm)

Catalog No.	Capacity (ml)
JG27901	1 ml
JG27902	10 ml

Note: Conical Flask with joint (Borosilicate Glass), Reagent Bottle (Borosilicate Glass) and Reagent Bottle (Neutral Glass) of 250/500ml available. Please specify at the time of quotation/order.

Note: Available with rubber cork in place of Nylon sleeve.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.