



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Tissue Flotation Bath

Product Image



Description

Tissue flotation bath is suitable for tissue culture, enzyme reactions fermentation, tissue section processing, etc. This is ideal for preparing distortion less and wrinkle free tissue specimens. **Construction:** It consists of a black powder coated aluminum tank housed in a powder coated outer case filled with high grade insulation provided with capillary type thermostat, indicator light ON/OFF switch and a lid. Temperature range ambient to 80°C +/- 1°C.

Particulars	Catalog No.
Tissue Floation Bath, 210 x 75mm Diameter x Depth	JA21410
Tissue Floation Bath, 250 x 125mm Diameter x Depth	JA21411

Note: Customized Baths (Water/Oil) with specific requirements like size, material, temperature range, etc are available on special demand

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.