



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Tubes, Culture (Disposable) Screw Thread

## Product Image

---



## Description

---

**Material:-**Borosilicate Glass

Tubes for Tissue Culture work.

**A type** Tubes, Culture (Disposable) Screw Thread, Flat Bottom with black phenolic cap with rubber liner

**B type** Tubes, Culture (Disposable) Screw Thread, Round Bottom with black phenolic cap with rubber liner

**C type** Tubes, Culture (Disposable) Screw Thread, Flat Bottom with white polypropylene screw cap

**D type** Tubes, Culture (Disposable) Screw Thread, Round Bottom with white polypropylene screw cap

**A type: Tubes, Culture (Disposable) Screw Thread, Flat Bottom with black phenolic cap with rubber liner**

<b>Catalog No.</b>	<b>O.D. x Length (mm)</b>	<b>Overflow Capacity (ml)</b>	<b>Case Qty.</b>
JG21201	13x100	8	100
JG21202	15x100	10	100
JG21203	16x100	12	100
JG21204	16x125	16	100
JG21205	16x150	20	100
JG21206	20x125	25	100
JG21207	20x150	30	100
JG21208	20x200	60	100

**B type: Tubes, Culture (Disposable) Screw Thread, Round Bottom with black phenolic cap with rubber liner**

<b>Catalog No.</b>	<b>O.D. x Length (mm)</b>	<b>Overflow Capacity (ml)</b>	<b>Case Qty.</b>
JG21301	13x100	8	100
JG21302	15x100	10	100
JG21303	16x100	12	100
JG21304	16x125	16	100
JG21305	16x150	20	100
JG21306	20x125	25	100
JG21307	20x150	30	100
JG21308	20x200	60	100

**C type: Tubes, Culture (Disposable) Screw Thread, Flat Bottom with white polypropylene screw cap**

<b>Catalog No.</b>	<b>O.D. x Length (mm)</b>	<b>Overflow Capacity (ml)</b>	<b>Case Qty.</b>
JG21401	13x100	8	100

JG21402	15x100	10	100
JG21403	16x100	12	100
JG21404	16x125	16	100
JG21405	16x150	20	100
JG21406	20x125	25	100
JG21407	20x150	30	100
JG21408	20x200	60	100

**D type: Tubes, Culture (Disposable) Screw Thread, Round Bottom with white polypropylene screw cap**

<b>Catalog No.</b>	<b>O.D. x Length (mm)</b>	<b>Overflow Capacity (ml)</b>	<b>Case Qty.</b>
JG21501	13x100	8	100
JG21502	15x100	10	100
JG21503	16x100	12	100
JG21504	16x125	16	100
JG21505	16x150	20	100
JG21506	20x125	25	100
JG21507	20x150	30	100
JG21508	20x200	60	100

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.