



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Universal Reagent Reservoir

Product Image



Description

Material: Polypropylene

This universal reagent reservoir is a graduated single reagent reservoir with maximum volume of 50 ml reservoir with 8 separate numbered 5 ml reagent channels and a reservoir with 12 separate numbered 5 ml reagent channel .Made of polypropylene these reagent reservoirs are repeatedly autoclavable and can be used with all multi channel pipettes

Catalog No.	Capacity (ml)
JW09901	50

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.