



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Universal Reagent Reservoir

## Product Image

---



## Description

---

**Material:** Polypropylene

This universal reagent reservoir is a graduated single reagent reservoir with maximum volume of 50 ml reservoir with 8 separate numbered 5 ml reagent channels and a reservoir with 12 separate numbered 5 ml reagent channel .Made of polypropylene these reagent reservoirs are repeatedly autoclavable and can be used with all multi channel pipettes

**Catalog No.**

JW09901

**Capacity (ml)**

50

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.