

Description

Catalog No. 100987

Urine Analyzer - Model 1240

TECHNICAL SPECIFICATIONS

High-luminance and white LED, features in long life and good stability.

Large LCD screen, high luminance, abundant contents display,

User-friendly interface.

Optional units: International unit, conventional unit and symbol system.

Three working mode: One-step/slow/fast testing mode,

Suitable for different user group.

Monitoring the whole test process, auto-character and audible prompt.

Be compatible with 8, 10 and 11-parameter test strip.

Standard RS232 interface and interface for data communication.

Built-in thermal printer.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.