



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Urinometer

Product Image



Description

Shot & wax poised, adjusted for use at 20°C. Specific Gravity Range: 1.000-1.060

Catalog No.	Model
JG27301	Urinometer (Black)
JG27302	Urinometer (Silver Coated)

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.