



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Weight Box

## Product Image



## Description

---

Material:- Polished brass

Precision polished brass masses (weights) in wooden case with hinged lid. It includes small plastic box which holds metal fractions. Complete with chrome plated brass forceps to lift masses. Each set consists of masses in 1-2-2-5 system in the standard range.

Catalog No.	Particulars
JZ03420	Weight Box, Analytical, Brass (1mg to 100g, Total 210g)
JZ03421	Weight Box, Physical, Brass (1mg to 100g, Total 210g)
JZ03422	Weight Box, Physical, Brass (1mg to 200g, Total 410g)
JZ03423	Weight Box, Physical, Brass (1mg to 500g, Total 1Kg)
JZ03424	Weight Box, Physical, Brass (1mg to 1000g, Total 2Kg)

A quality product from JAPSON ideal for use in School Laboratories in developing countries of Asia, Africa and South America.

Standard specifications for school laboratory use. Ideal mix of Quality and Economy. Other configurations and superior quality available for specific use or for developed countries.

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.