



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# pH Meter Scale, Bench Type

## Product Image



## Description

---

Highly Stable & Accurate  
Cell Constant Adjustment Facility  
3.5 Digit Red LED Display

### Technical Specifications :

General

Power: 220 V AC  $\pm$  10%, 50Hz

### Parameters:

pH Range: 0 - 14 pH

Resolution: 0.01 pH

Accuracy:  $\pm$  0.01 pH  $\pm$  1 Digit

Temperature Compensation Manual: 0 to 100°C

Temperature

Range: 0 - 100°C

Resolution: 0.1°C

Accuracy:  $\pm$  0.2%  $\pm$  1 Digit

Sensor RTD: (PT - 100)

Conductivity

Range: 0 to 20 mS/cm in two Ranges

Accuracy:  $\pm$  0.5%  $\pm$  1 Digit

Resolution: 1 $\mu$ S/cm

### Accessories :

pH combination electrode

Temperature probe

Operation manual

Dust cover

pH stand

7 & 4 pH bottles

Beaker

Catalog No.	Particulars
-------------	-------------

JZ02400	pH Meter, Bench Type
---------	----------------------

A quality product from JAPSON ideal for use in School Laboratories in developing countries of Asia, Africa and South America.

Standard specifications for school laboratory use. Ideal mix of Quality and Economy.

Other configurations and superior quality available for specific use or for developed

countries.

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.